

# Oticon Connect Clip Manual

## ConnectClip Quick Start Guide

### PAIRING CONNECTCLIP WITH HEARING AIDS

1. Ensure ConnectClip is powered off. Press and hold multi-button for approximately 6 seconds until power indicator OFF ●
2. Turn both hearing aids OFF. (Open battery drawer)
3. Turn both hearing aids ON. (Close battery drawer)
4. The hearing aids will now be ready for first time connection with the ConnectClip. You now have 3 minutes to complete the pairing.
5. Turn ON ConnectClip. Press and hold multi-button for approximately 3 seconds until power indicator shows STEADY GREEN ●
6. ConnectClip will now search for



### PAIRING CONNECTCLIP WITH MOBILE PHONE

1. Pair the ConnectClip with your hearing aids before pairing ConnectClip with your mobile phone.
2. Place ConnectClip in Pairing mode. Press and hold the multi-button AND volume up button for about 6 seconds (A) until status BLINKS BLUE ●
3. Enable pairing on mobile phone. Turn on Bluetooth. Search for and connect to/pair with ConnectClip. This is normally in the Settings/Bluetooth section of your phone.
4. Enter your pin code, if your phone



**Oticon Connect Clip Manual** is an essential guide for users of Oticon hearing aids who want to enhance their auditory experience. The Connect Clip is a versatile accessory that allows users to connect their hearing aids to smartphones, tablets, and other Bluetooth-enabled devices. In this article, we will explore the features, setup process, usage tips, and troubleshooting for the Oticon Connect Clip.

## Understanding the Oticon Connect Clip

The Oticon Connect Clip is designed to improve communication and accessibility for individuals with hearing impairment. This device acts as a bridge between your hearing aids and various audio sources, providing a seamless listening experience.

## Key Features of the Oticon Connect Clip

The Oticon Connect Clip comes with several notable features:

- **Bluetooth Connectivity:** Connects to various devices, including smartphones, tablets, and computers.
- **Remote Microphone:** Allows users to hear conversations clearly in noisy environments.
- **Streaming Audio:** Enables streaming of phone calls, music, and other audio directly to the hearing aids.
- **Hands-Free Functionality:** Acts as a hands-free device for phone calls, making communication easier.
- **Multi-Device Pairing:** Can be connected to multiple devices simultaneously, allowing for easy switching between them.

# Setting Up Your Oticon Connect Clip

Setting up your Oticon Connect Clip is a straightforward process. Follow these steps to get started:

## 1. Charging the Connect Clip

Before using the Connect Clip, ensure that it is fully charged:

- Use the provided USB charging cable to connect the device to a power source.
- A full charge typically takes about 2 hours and provides several hours of usage.

## 2. Pairing the Connect Clip with Your Hearing Aids

To pair the Connect Clip with your Oticon hearing aids:

1. Turn on the Hearing Aids: Ensure that your hearing aids are powered on and within range of the Connect Clip.
2. Activate Pairing Mode: Press and hold the Connect Clip button until the LED indicator flashes blue.
3. Use the Oticon ON App: If you have the Oticon ON app, open it and follow the on-screen instructions to complete pairing.
4. Confirm Connection: Once paired, the LED indicator will change to a steady blue light.

## 3. Pairing with Bluetooth Devices

To pair the Connect Clip with Bluetooth-enabled devices:

1. Enable Bluetooth on the Device: Go to the settings of your smartphone or tablet and enable Bluetooth.
2. Activate Pairing Mode on Connect Clip: Press and hold the Connect Clip button until the LED flashes blue.
3. Select the Connect Clip: On your device, find and select "Oticon Connect Clip" from the list of available devices.
4. Confirm Connection: Once paired, the LED will turn to a steady blue light.

## Using the Oticon Connect Clip

Once you have set up the Connect Clip, you can begin to take advantage of its features.

# **1. Streaming Audio**

To stream audio such as phone calls or music:

- Ensure the Connect Clip is paired with your Bluetooth device.
- Start playing audio on your device, and it will stream directly to your hearing aids.

# **2. Using the Remote Microphone Feature**

The remote microphone feature enhances the ability to hear conversations in noisy environments:

- Place the Connect Clip on the table or give it to the person speaking.
- Press the button on the Connect Clip to activate the microphone.
- Audio from the person speaking will be transmitted directly to your hearing aids, minimizing background noise.

# **3. Making Phone Calls**

To make phone calls using the Connect Clip:

- Use your smartphone to dial the number.
- The audio will stream through the Connect Clip, allowing for hands-free conversation.
- Speak into the Connect Clip's microphone for clear communication.

# **Maintaining Your Oticon Connect Clip**

Proper maintenance of your Connect Clip will ensure its longevity and optimal performance:

## **1. Regular Cleaning**

- Wipe the device with a soft, dry cloth to remove dirt and moisture.
- Avoid using cleaning agents or solvents that can damage the device.

## **2. Battery Care**

- Avoid letting the battery completely drain; charge it regularly to maintain battery health.
- Store the Connect Clip in a cool, dry place when not in use.

# Troubleshooting Common Issues

Despite its ease of use, users may encounter issues with the Oticon Connect Clip. Here are some common problems and solutions:

## 1. Connect Clip Won't Pair with Hearing Aids

- Solution: Ensure that both the hearing aids and Connect Clip are turned on and within range. Restart both devices and try pairing again.

## 2. Poor Audio Quality During Streaming

- Solution: Check for obstacles between the Connect Clip and your audio device. Ensure that both devices are fully charged.

## 3. Connect Clip Not Charging

- Solution: Inspect the charging cable and port for any damage. Try a different power source or cable if necessary.

## Conclusion

The **Oticon Connect Clip Manual** serves as a crucial resource for users looking to maximize their experience with Oticon hearing aids. By understanding the features, following the setup process, and utilizing the various functions of the Connect Clip, users can enjoy enhanced communication and audio streaming. Regular maintenance and troubleshooting can further ensure that the device remains in good working condition, providing long-term support for individuals with hearing impairments. With the Oticon Connect Clip, users can experience a more connected and enriching auditory world.

## Frequently Asked Questions

### What is the purpose of the Oticon ConnectClip?

The Oticon ConnectClip serves as a wireless accessory that allows users to connect their hearing aids to various devices such as smartphones, tablets, and other audio sources, enhancing sound quality and convenience.

## Where can I find the Oticon ConnectClip manual?

The Oticon ConnectClip manual can be found on the official Oticon website under the support section, or by contacting your hearing care professional for a printed copy.

## How do I pair my Oticon ConnectClip with my smartphone?

To pair your Oticon ConnectClip with your smartphone, turn on Bluetooth on your phone, press and hold the ConnectClip's button until the pairing mode is activated, then select 'ConnectClip' from the list of available devices on your smartphone.

## What should I do if my Oticon ConnectClip is not charging?

If your Oticon ConnectClip is not charging, check the charging cable and power source for any issues, ensure the ConnectClip is properly seated in the charging dock, and clean the charging contacts to remove any dirt or debris.

## Can I use the Oticon ConnectClip for phone calls?

Yes, the Oticon ConnectClip is designed to facilitate phone calls directly through your hearing aids, providing clearer sound and hands-free functionality.

## What are the main features of the Oticon ConnectClip?

Main features of the Oticon ConnectClip include Bluetooth connectivity, the ability to stream audio from devices, hands-free calling, and access to remote microphone functionality for improved listening in noisy environments.

## How do I reset my Oticon ConnectClip?

To reset your Oticon ConnectClip, press and hold the button for about 10 seconds until the LED indicator blinks red and blue, then release the button. This will restore factory settings and clear any paired devices.

Find other PDF article:

<https://soc.up.edu.ph/20-pitch/files?trackid=ElS74-2473&title=everyday-ways-to-lose-weight.pdf>

# Oticon Connect Clip Manual

□ □ □ □ □ □ □ □ □ □ □ □ □ □ □ □ □ □ - □ □

Oticon 100 Demant ...

助聴器の種類と特徴 - 助聴器

Oticon/Bernafron 助聴器の種類と特徴 - 助聴器

助聴器の種類と特徴 - 助聴器

助聴器の種類と特徴 PHONAK 助聴器 SONOVA 助聴器 Oticon 助聴器 Demant 助聴器 Resound 助聴器 Starkey 助聴器 ...

助聴器の種類と特徴 - 助聴器

助聴器の種類と特徴 4 助聴器の種類と特徴 1 助聴器の種類と特徴 ...

助聴器の種類と特徴? - 助聴器

助聴器の種類と特徴 PHONAK 助聴器 SONOVA 助聴器 Oticon 助聴器 Demant 助聴器 Resound 助聴器 Starkey ...

助聴器の種類と特徴 - 助聴器

助聴器の種類と特徴 “助聴器” 助聴器の種類と特徴 ...

助聴器の種類と特徴 - 助聴器

Oct 29, 2020 · 助聴器の種類と特徴 Bone Anchored ...

助聴器の種類と特徴 - 助聴器

助聴器の種類と特徴 ...

助聴器の種類と特徴 - 助聴器

助聴器の種類と特徴 Oticon 助聴器 100 助聴器の種類と特徴 Demant 助聴器 ...

助聴器の種類と特徴 - 助聴器

助聴器の種類と特徴 Oticon/Bernafron 助聴器の種類と特徴 - 助聴器

助聴器の種類と特徴 - 助聴器

助聴器の種類と特徴 PHONAK 助聴器 SONOVA 助聴器 Oticon 助聴器 Demant 助聴器 Resound 助聴器 Starkey 助聴器 WS Audiology 助聴器 Signia 助聴器 WIDEX 助聴器

助聴器の種類と特徴 - 助聴器

助聴器の種類と特徴 4 助聴器の種類と特徴 1 助聴器の種類と特徴 ...

助聴器の種類と特徴? - 助聴器

助聴器の種類と特徴 PHONAK 助聴器 SONOVA 助聴器 Oticon 助聴器 Demant 助聴器 Resound 助聴器 Starkey 助聴器 Signia 助聴器 WS Audiology 助聴器 WIDEX 助聴器 WS Audiology 助聴器 2019 WIDEX ...

助聴器の種類と特徴 - 助聴器

助聴器の種類と特徴 “助聴器” 助聴器の種類と特徴 ...

骨導式人工耳蜗 - 骨

Oct 29, 2020 · 骨導式人工耳蜗 (BAHA) 是一种将声音信号通过头骨传导到内耳的设备。Bone Anchored Hearing Aid (BAHA) 是一种将声音信号通过头骨传导到内耳的设备。

骨導式人工耳蜗 - 骨

骨導式人工耳蜗 (BAHA) 是一种将声音信号通过头骨传导到内耳的设备。骨導式人工耳蜗 (BAHA) 是一种将声音信号通过头骨传导到内耳的设备。 ...

Unlock the full potential of your hearing experience with the Oticon Connect Clip manual. Learn how to maximize its features and enhance your connectivity today!

[Back to Home](#)